

Technical Data Sheet

Product

Anti-Hu IgD Purified

Cat. Number/Size

11-142-C025

0.025 mg

11-142-C100

0.1 mg

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen

IgD

Clone

IA6-2

Format

Purified

Reactivity

Human

Application

FC (QC tested), IHC(F)

Application details

Flow cytometry: Recommended dilution: 2-4 µg/ml.

Isotype

Mouse IgG2a kappa

Specificity

The mouse monoclonal antibody IA6-2 recognizes human immunoglobulin D.

Other names

Immunoglobulin D

Immunogen

Human IgD

Concentration

1 mg/ml

Preparation

Purified by protein-A affinity chromatography.

Formulation

Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage and handling

Store at 2-8°C. Do not freeze.

Images and References

www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.