

Technical Data Sheet

Product	Anti-Neurofilament M Purified	
Cat. Number/Size	11-294-C025	0.025 mg
	11-294-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	Neurofilament M
Clone	NF-09
Format	Purified
Reactivity	Mammalian
Application	IHC(P), WB, ICC
Isotype	Mouse IgG2a
Specificity	The antibody NF-09 reacts with both phosphorylated and non-phosphorylated form of medium neurofilament protein (160 kDa intracellular antigen) of various species.
Other names	NFM, NEFM, Neurofilament medium chain, NEF3
Immunogen	Pellet of porcine brain cold stable proteins after depolymerization of microtubules.
Entrez Gene ID	4741
Gene name	NEFM
NCBI Full Gene Name	neurofilament medium
UniProt ID	P07197

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.