

Technical Data Sheet

Product	Anti-Hu STAT1 (pS727) Purified	
Cat. Number/Size	11-372-C025	0.025 mg
	11-372-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	STAT1 (pS727)
Clone	PSM1
Format	Purified
Reactivity	Human
Application	IP, WB
Isotype	Mouse IgG1
Specificity	The antibody PSM1 recognizes transcriptional factor STAT1 (91 kDa) activated by phosphorylation at Ser727 (intracellular antigen).
Other names	CANDF7, IMD31, ISGF3
Immunogen	STAT1 peptide sequence 721-733 (Ser727 phosphorylated).
Entrez Gene ID	6772
Gene name	STAT1
NCBI Full Gene Name	signal transducer and activator of transcription 1
UniProt ID	P42224

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.