



	Technical Data Sheet	
Product	Anti-Hu CD230 / Prion Purified	
Cat. Number/Size	11-593-C025	0.025 mg
	11-593-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD230 / Prion	
Clone	EM-20	
Format	Purified	
Reactivity	Human	
Application	WB	
Application details	Western blotting: Recommended dilution: 0.5 µg/ml, non-reducing conditions are essential.	
Isotype	Mouse IgG2a	
Specificity	The mouse monoclonal antibody EM-20 recognizes human prion protein (PrP). Diglycosylated form of PrP has ~ 40 kDa, monoglycosylated form ~ 30 kDa, and nonglycosylated form ~ 19-21 kDa. This antibody is suitable for discrimination between normal cellular prion protein (PrPc) and its conformationally changed form (PrPSc) prion protein.	
Other names	Major prion protein, PrP, PrP27-30, PrP33-35C, ASCR, CD230, PRNP, CJD, KURU, PRIP, PRPC	
Immunogen	Recombinant human prion protein	
Entrez Gene ID	5621	
Gene name	PRNP	
NCBI Full Gene Name	prion protein	
UniProt ID	P04156	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.