

Technical Data Sheet

Product	Anti-SLP76 Purified	
Cat. Number/Size	11-402-C025	0.025 mg
	11-402-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	SLP76
Clone	SLP-76/03
Format	Purified
Reactivity	Pig, Human, Mouse
Application	WB, IHC-P
Application details	Immunohistochemistry (paraffin sections): Recommended dilution: 10 µg/ml; positive tissue: thymus
Isotype	Mouse IgG2b
Specificity	The monoclonal antibody reacts with SLP76, a 76kDa cytosolic adaptor protein that is involved in signaling of various hematopoietic cells, such as T cells, mast cells or neutrophils; in B cells, however, it is replaced by SLP65.
Other names	SLP-76, LCP2
Immunogen	Bacterially expressed fusion protein representing amino acids 216-434 of human SLP76 with histidine tag
Entrez Gene ID	3937
Gene name	LCP2
NCBI Full Gene Name	lymphocyte cytosolic protein 2
UniProt ID	Q13094

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.