

Technical Data Sheet

Product	Anti-Hu CD3 zeta (pY72) Purified	
Cat. Number/Size	11-712-C025	0.025 mg
	11-712-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD3 zeta (pY72)
Clone	EM-26
Format	Purified
Reactivity	Mouse, Human
Application	FC, WB, ICC
Isotype	Mouse IgG2b
Specificity	The mouse monoclonal antibody EM-26 recognizes phosphorylated tyrosine 72 of CD3 zeta chain (CD247), which is a component of TCR/CD3 complex expressed on T cells.
Other names	CD3 zeta chain, CD247, T3Z, TCRz, IMD25
Immunogen	A phospho specific peptide corresponding to the amino acids surrounding tyrosine 72 of mouse CD3 zeta linked to KLH

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on vial label.

Images and References www.exbio.cz

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic use. In vivo diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s. © 2019 EXBIO Praha, a.s. All rights reserved.