

Technical Data Sheet

Product	Anti-Hu CD117 APC	
Cat. Number/Size	1A-586-T025	25 tests
	1A-586-T100	100 tests
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD117	
Clone	104D2	
Format	APC	
Reactivity	Non-human primates, Human, Cow	
Application	FC (QC tested)	
Application details	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.	
Excitation laser	red (633 nm)	
Isotype	Mouse IgG1 kappa	
Specificity	The mouse monoclonal antibody 104D2 detects extracellular part of CD117 / c-Kit protooncogen.	
Other names	c-Kit, PBT, SCFR	
Workshop	HLDA VI: WS Code C-30	
Immunogen	MOLM-1 megakaryocytic cells	
Entrez Gene ID	3815	
Gene name	KIT	
NCBI Full Gene Name	KIT proto-oncogene, receptor tyrosine kinase	
UniProt ID	P10721	
Preparation	The purified antibody is conjugated with allophycocyanin (APC) under optimum conditions. The conjugate is purified by size-exclusion chromatography.	
Formulation	Stabilizing phosphate buffered saline (PBS) solution containing 15 mM sodium azide	
Storage and handling	Store in the dark at 2-8°C. Do not freeze. Avoid prolonged exposure to light. Do not use after expiration date stamped on vial label.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.