



Technical Data Sheet

Product	Anti-Hu Tomm34 Purified	
Cat. Number/Size	11-891-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	Tomm34	
Clone	Tomm34-4.1	
Format	Purified	
Reactivity	Human	
Application	WB, ICC, IP	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody Tomm34-4.1 recognizes an epitope within TPR2 domain of human Tomm34 protein.	
Other names	URCC3, HTOM34P	
Immunogen	Recombinant human Tomm34	
Entrez Gene ID	10953	
Gene name	TOMM34	
NCBI Full Gene Name	Translocase of outer mitochondrial membrane 34	
UniProt ID	Q15785	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.