



	Technical Data Sheet	
Product	Anti-HSP90 alpha/beta Purified	
Cat. Number/Size	11-533-C025	0.025 mg
	11-533-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	HSP90 alpha/beta	
Clone	MBH90AB	
Format	Purified	
Reactivity	Cow, Mouse, Human	
Application	WB, IHC-P	
Application details	Western blotting: Recommended dilution: 1 µg/ml.	
Isotype	Mouse IgG1	
Specificity	The antibody MBH90AB recognizes the epitope EEEVE within N-terminal part of ubiquitously expressed Hsp90 alpha and Hsp90 beta intracellular proteins with calculated Mw of 84.7 kDa and 83.3 kDa, respectively, however, migrating as 90 kDa bands under reducing SDS-PAGE conditions.	
Other names	HSP90A, HSP90B, LAP-2, LAP2, EL52	
Immunogen	Peptide corresponding to the sequence EEVHHGEEVEEC within N-terminal part of human Hsp90.	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.