

Technical Data Sheet

Product	<b>Anti-DLL4 Purified</b>	
Cat. Number/Size	<b>11-717-C025</b>	<b>0.025 mg</b>
	<b>11-717-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	DLL4
Clone	MHD4-46
Format	Purified
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody MHD4-46 recognizes the extracellular domain of DLL4 (Delta-like ligand 4), a type I transmembrane protein which plays an important role in vascular development.
Other names	Delta like ligand 4, AOS6, canonical Notch ligand 4
Immunogen	recombinant soluble human DLL4
Entrez Gene ID	54567
Gene name	DLL4
NCBI Full Gene Name	delta like canonical Notch ligand 4
UniProt ID	Q9NR61

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.