



	Technical Data Sheet	
Product	Anti-Hu CD300e Purified	
Cat. Number/Size	11-728-C025	0.025 mg
	11-728-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD300e	
Clone	UP-H2	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), IP, WB	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.	
Isotype	Mouse IgG1 kappa	
Specificity	The mouse monoclonal antibody UP-H2 recognizes an extracellular epitope on CD300e / IREM-2, a 32 kDa glycoprotein expressed by mature monocytes and peripheral blood myeloid dendritic cells.	
Other names	CLM2, CMRF35-A5, LMIR6, IREM-2, PIgR2	
Immunogen	CD300e-HA-transfected cells	
Entrez Gene ID	342510	
Gene name	CD300E	
NCBI Full Gene Name	CD300e molecule	
UniProt ID	Q496F6	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.