



	Technical Data Sheet	
Product	Anti-ARAP1 Purified	
Cat. Number/Size	11-516-C025	0.025 mg
	11-516-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	ARAP1	
Clone	ARAP1-2	
Format	Purified	
Reactivity	Human	
Application	IP, WB, ICC	
Application details	Western blotting: Recommended dilution: 2 µg/ml	
Isotype	Mouse IgG1	
Specificity	The antibody ARAP1-2 reacts with C-terminal part of ARAP1 (intracellular epitope), a 160 kDa adaptor protein.	
Other names	CED2, CENTD2	
Immunogen	His6-tagged protein corresponding to amino acids 1190-1450 of human ARAP1	
Entrez Gene ID	116985	
Gene name	ARAP1	
NCBI Full Gene Name	ArfGAP with RhoGAP domain, ankyrin repeat and PH d	
UniProt ID	Q96P48	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.