



	Technical Data Sheet	
Product	Anti-PCNA Purified	
Cat. Number/Size	11-182-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	PCNA	
Clone	PC10	
Format	Purified	
Reactivity	Mouse, Human, Non-human primates, Drosophila, Chicken, Rat	
Application	FC (QC tested), IP, WB, IHC-P, IHC-F, ICC	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.	
Isotype	Mouse IgG2a kappa	
Specificity	The mouse monoclonal antibody PC10 (also known as 3F81) recognizes PCNA, a 36 kDa conserved nuclear protein serving as a cofactor for DNA synthesis.	
Other names	ATLD2	
Immunogen	recombinant rat PCNA	
Entrez Gene ID	5111	
Gene name	PCNA	
NCBI Full Gene Name	proliferating cell nuclear antigen	
UniProt ID	P12004	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.