



	Technical Data Sheet	
Product	Anti-GRAP2 APC	
Cat. Number/Size	1A-639-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	GRAP2	
Clone	UW40	
Format	APC	
Reactivity	Human	
Application	FC-IC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.	
Excitation laser	red (633 nm)	
Isotype	Mouse IgG2a	
Specificity	The mouse monoclonal antibody UW40 recognizes GRAP2/GADS, a 41 kDa cytoplasmic adaptor protein that plays a pivotal role during the early events of signal transduction in T cells.	
Other names	GADS, GRID, GRAP-2, GrbX	
Immunogen	GST-fusion human GRAP2/GADS protein	
Entrez Gene ID	9402	
Gene name	GRAP2	
NCBI Full Gene Name	GRB2 related adaptor protein 2	
UniProt ID	O75791	
Concentration	0.1 mg/ml	
Preparation	Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.	
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.