

Technical Data Sheet

Product

Anti-GST tag Purified

Cat. Number/Size

11-477-C025

0.025 mg

11-477-C100

0.1 mg

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	GST tag
Clone	S-tag-05
Format	Purified
Reactivity	Tagged fusion proteins in all species
Application	WB, IP
Application details	Western blotting: Recommended dilution: 1 µg/ml; positive control: GST fusion protein, reducing conditions.
Isotype	Mouse IgG2b
Specificity	The antibody S-tag-05 is a mouse monoclonal antibody with high affinity to the glutathione-S-transferase (GST) from <i>Schistosoma japonicum</i> . This purified antibody is suitable for detecting fusion proteins containing the GST-Tag sequence expressed in <i>E. coli</i> , yeast, mammalian, and <i>in vitro</i> transcription/translation systems.
Immunogen	GST fusion protein

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.