



	Technical Data Sheet	
Product	Anti-Ki-67 Purified	
Cat. Number/Size	11-155-C025	0.025 mg
	11-155-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	Ki-67	
Clone	Ki-67	
Format	Purified	
Reactivity	Cow, Human	
Application	FC (QC tested), WB, IHC, ICC	
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml. Intracellular staining. Immunocytochemistry: Paraformaldehyde fixation; recommended antibody concentration 1 µg/ml.	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody Ki-67 recognizes Ki-67 antigen, a non-histone nuclear protein expressed exclusively in proliferating cells.	
Other names	KIA, MIB, Marker of proliferation Ki67, MKI67	
Immunogen	Nuclei of the Hodgkin lymphoma cell line L428	
Entrez Gene ID	4288	
Gene name	MKI67	
NCBI Full Gene Name	marker of proliferation Ki-67	
UniProt ID	P46013	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.