



	Technical Data Sheet	
Product	Anti-Hu CD88 PE	
Cat. Number/Size	1P-736-T025	25 tests
	1P-736-T100	100 tests
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD88	
Clone	S5/1	
Format	PE	
Reactivity	Rabbit, Ferret, Cow, Human	
Application	FC (QC tested)	
Application details	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.	
Excitation laser	blue (488 nm)	
Isotype	Mouse IgG2a kappa	
Specificity	The mouse monoclonal antibody S5/1 recognizes an extracellular epitope of CD88 protein, a 43 kDa receptor of C5a component of the complement cascade.	
Other names	C5A, C5R1, C5AR1	
Workshop	HCDM listed	
Immunogen	Recombinant N-terminal peptide (Asp15-Asp27) of human C5aR	
Entrez Gene ID	728	
Gene name	C5AR1	
NCBI Full Gene Name	complement C5a receptor 1	
UniProt ID	P21730	
Preparation	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.	
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.