

Technical Data Sheet

Product	Anti-EBNA-1 (EBV) Purified	
Cat. Number/Size	11-303-C025	0.025 mg
	11-303-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	EBNA-1 (EBV)
Clone	E1-2.5
Format	Purified
Reactivity	Epstein-Barr Virus
Application	ICC, ELISA
Isotype	Mouse IgG2b
Specificity	The mouse monoclonal antibody E1-2.5 recognizes repetitive Gly-Ala region of the Epstein-Barr virus (EBV) nuclear antigen 1 (EBNA-1). EBNA-1 is the only viral protein consistently expressed in all malignancies associated with EBV.
Entrez Gene ID	3783709
Gene name	EBNA-1
NCBI Full Gene Name	nuclear antigen EBNA-1
UniProt ID	P03211

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.