



	Technical Data Sheet	
Product	<b>Anti-Hu CD50 Purified</b>	
Cat. Number/Size	<b>11-266-C025</b>	<b>0.025 mg</b>
	<b>11-266-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	
Antigen	CD50	
Clone	MEM-171	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), IP	
Application details	Flow cytometry: Recommended dilution: 4 µg/ml.	
Isotype	Mouse IgG1	
Specificity	The antibody MEM-171 recognizes an extracellular epitope in the D2 domain of CD50 (ICAM-3), a 120-130 kDa type I membrane protein (immunoglobulin supergene family) expressed on leukocytes, endothelial cells and Langerhans cells; it is negative on platelets and erythrocytes.	
Other names	ICAM-3, ICAM-R	
Workshop	HLDA VI: WS Code BP 614; HLDA VI: WS Code NL N-L022	
Immunogen	Human granulocytes	
Entrez Gene ID	3385	
Gene name	ICAM3	
NCBI Full Gene Name	intercellular adhesion molecule 3	
UniProt ID	P32942	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	<a href="http://www.exbio.cz">www.exbio.cz</a>	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.