

Technical Data Sheet

Product

Anti-H-RAS Purified

Cat. Number/Size

11-527-C025

0.025 mg

11-527-C100

0.1 mg

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	H-RAS
Clone	H-RAS-03
Format	Purified
Reactivity	Human
Application	WB
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody H-RAS-03 reacts with human H-Ras, an ubiquitously expressed 21 kDa intracellular protein. Although reactivity with other species has not been determined, it is probable as the epitope is highly conserved among animals.
Other names	RAS-H, HRAS, p21Ras
Immunogen	Peptide corresponding to amino acids DIHQYREQIKRVKDSDDC of human H-Ras protein
Entrez Gene ID	3265
Gene name	HRAS
NCBI Full Gene Name	HRas proto-oncogene
UniProt ID	P01112

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.