



	Technical Data Sheet	
Product	Anti-Hu CD62L Purified	
Cat. Number/Size	11-630-C025	0.025 mg
	11-630-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD62L	
Clone	DREG56	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), IP, WB, IHC-F	
Application details	Flow cytometry: Recommended dilution: 2 µg/ml. Western blotting: Non-reducing conditions.	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody DREG56 recognizes an extracellular epitope of CD62L / L-selectin, a 65-76 kDa cell surface protein, expressed by neutrophils, monocytes, and subsets of T, B, and NK cells, that interacts with specific carbohydrates exposed on activated endothelial cells.	
Other names	L-selectin, LSEL, SELL, LAM-1, LECAM1, LEU8, LNHR, LYAM1, PLNHR	
Workshop	HLDA V: WS Code S056	
Immunogen	PMA-activated human peripheral blood leukocytes	
Entrez Gene ID	6402	
Gene name	SELL	
NCBI Full Gene Name	selectin L	
UniProt ID	P14151	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.