



	Technical Data Sheet	
Product	Anti-Hu CD326 Biotin	
Cat. Number/Size	1B-582-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD326	
Clone	323/A3	
Format	Biotin	
Reactivity	Human	
Negative species	Rat	
Application	FC (QC tested), IHC-P, WB, ICC, IP	
Application details	Flow cytometry: Recommended dilution: 0.5-4 µg/ml.	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody 323/A3 recognizes an extracellular epitope of CD326 / EpCAM, a marker of epithelial lineages, that is over-expressed in many carcinomas.	
Other names	EPCAM, GA733-2, EGP314, KSA, KS 1/4 antigen, Trop-1, M4S1, DIAR5, MIC18, TROP1	
Immunogen	Human breast cancer MCF-7 cells	
Entrez Gene ID	4072	
Gene name	EPCAM	
NCBI Full Gene Name	epithelial cell adhesion molecule	
UniProt ID	P16422	
Concentration	1 mg/ml	
Preparation	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.