



	Technical Data Sheet	
Product	Anti-p63 Purified	
Cat. Number/Size	11-888-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	p63	
Clone	PANp63-6.1	
Format	Purified	
Reactivity	Mouse, Rat, Human	
Application	WB, IHC-P, ICC, IP	
Application details	Western blotting: Recommended dilution: 1 µg/ml	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody PANp63-6.1 recognizes both TAp63, and deltaNp63 form of p63. The target epitope PSHLIR is located within amino acids 261-266 of TAp63, and 167-172 of deltaNp63.	
Other names	TAp63 + deltaNp63, AIS, KET, LMS, NBP, RHS, EEC3, OFC8, SHFM4	
Immunogen	Recombinant full length human p63	
Entrez Gene ID	8626	
Gene name	TP63	
NCBI Full Gene Name	tumor protein p63	
UniProt ID	Q9H3D4	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.