



	Technical Data Sheet	
Product	Anti-HLA-G PE	
Cat. Number/Size	1P-431-C025	0.025 mg
	1P-431-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	HLA-G	
Clone	5A6G7	
Format	PE	
Reactivity	Human	
Application	FC-IC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.	
Excitation laser	blue (488 nm)	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody 5A6G7 was generated to a peptide corresponding to C-intron 4-encoded sequence. This antibody does not crossreact with the full-length HLA-G1 isoform and thus allows to distinguish between secreted HLA-G5 and HLA-G6 isoforms from shedded HLA-G1.	
Immunogen	C-terminal amino acid sequence (22-mer) of soluble HLA-G5 and HLA-G6 proteins coupled to ovalbumin.	
Entrez Gene ID	3135	
Gene name	HLA-G	
NCBI Full Gene Name	major histocompatibility complex, class I, G	
UniProt ID	P17693	
Concentration	0.1 mg/ml	
Preparation	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.	
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.