



	Technical Data Sheet	
Product	Anti-Hu NKp80 Purified	
Cat. Number/Size	11-917-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	NKp80	
Clone	5D12	
Format	Purified	
Reactivity	Human, Monkey	
Negative species	Mouse, Rat	
Application	FC (QC tested), IP	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.	
Isotype	Mouse IgG1 kappa	
Specificity	The mouse monoclonal antibody 5D12 recognizes an extracellular epitope of human NKp80 (CLEC5C), a C-type lectin family member, expressed on NK cells and subsets of T cells.	
Other names	KLRF1, CLEC5C	
Immunogen	recombinant human NKp80 extracellular domain	
Entrez Gene ID	51348	
Gene name	KLRF1	
NCBI Full Gene Name	killer cell lectin like receptor F1	
UniProt ID	Q9NZS2	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.