

Technical Data Sheet

Product

**Anti-Hu CD10 Purified**

Cat. Number/Size

**11-209-C025**

**0.025 mg**

**11-209-C100**

**0.1 mg**

**For Research Use Only.**

**Not for use in diagnostic or therapeutic procedures.**

---

Antigen	CD10
Clone	MEM-78
Format	Purified
Reactivity	Human
Application	FC (QC tested), IP
Application details	Flow cytometry: Recommended dilution: 1 µg/ml.
Isotype	Mouse IgG1
Specificity	The antibody MEM-78 reacts with an extracellular epitope CD10 antigen (CALLA - Common acute lymphatic leukemia antigen), a 100 kDa type II integral membrane protein.
Other names	CALLA, Neprilysin, Neutral endopeptidase, Enkephalinase, Atriopeptidase, MME
Workshop	HLDA IV: WS Code B 506; HLDA V: WS Code B CD10.4
Immunogen	NALM-6 human pre-B cell line
Entrez Gene ID	4311
Gene name	MME
NCBI Full Gene Name	membrane metalloendopeptidase
UniProt ID	P08473

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.