

Technical Data Sheet

Product

Anti-Hu CD157 APC

Cat. Number/Size

1A-136-T025

25 tests

1A-136-T100

100 tests

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	CD157
Clone	SY11B5
Format	APC
Reactivity	Non-human primates, Human
Application	FC (QC tested)
Application details	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Excitation laser	red (633 nm)
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody SY11B5 recognizes CD157, an approximately 45 kDa GPI-anchored extracellular protein expressed mainly on monocytes, macrophages, granulocytes and bone marrow stromal cells.
Other names	BST1, cADPr hydrolase 2
Workshop	HCDM listed
Immunogen	Human CD157
Entrez Gene ID	683
Gene name	BST1
NCBI Full Gene Name	bone marrow stromal cell antigen 1
UniProt ID	Q10588

Preparation	Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.