

Technical Data Sheet

Product	Anti-TNF alpha PE	
Cat. Number/Size	1P-688-T025	25 tests
	1P-688-T100	100 tests
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	TNF alpha
Clone	MAb11
Format	PE
Reactivity	Pig, Human, Non-human primates
Application	FC (QC tested)
Application details	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Excitation laser	blue (488 nm)
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody MAb11 recognizes human 17-26 kDa cytokine TNF-alpha (tumor necrosis factor alpha).
Other names	TNFA, DIF, TNLG1F
Immunogen	Recombinant human TNF-alpha
Entrez Gene ID	7124
Gene name	TNF
NCBI Full Gene Name	tumor necrosis factor
UniProt ID	P01375

Preparation	The purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. The conjugate is purified by size-exclusion chromatography.
Formulation	Stabilizing phosphate buffered saline (PBS) solution containing 15 mM sodium azide
Storage and handling	Store in the dark at 2-8°C. Do not freeze. Avoid prolonged exposure to light. Do not use after expiration date stamped on vial label.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.