



	Technical Data Sheet	
Product	Anti-TFG Purified	
Cat. Number/Size	11-709-C025	0.025 mg
	11-709-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	TFG	
Clone	TFG-03	
Format	Purified	
Reactivity	Human	
Application	IP, WB, ICC	
Isotype	Mouse IgG2b	
Specificity	The mouse monoclonal antibody TFG-03 recognizes TFG, an approximately 50 kDa intracellular protein with regulatory functions.	
Other names	TF6, TRKT3, SPG57	
Immunogen	Human recombinant protein Trx-his-NTFG	
Entrez Gene ID	10342	
Gene name	HMSNP	
NCBI Full Gene Name	trafficking from ER to golgi regulator	
UniProt ID	Q92734	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.