

Technical Data Sheet

Product	Anti-Hu IgG (Fab) Purified	
Cat. Number/Size	11-319-C025	0.025 mg
	11-319-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	IgG (Fab)	
Clone	4A11	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), WB, ELISA	
Isotype	Mouse IgG2b	
Specificity	The mouse monoclonal antibody 4A11 reacts with Fab fragment of human IgG.	
Other names	Immunoglobulin G Fab fragment	
Immunogen	Purified human IgG.	
Concentration	1 mg/ml	
Preparation	Purified by precipitation and chromatography	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.