

Technical Data Sheet

Product **Anti-Hu TCR Valpha24-Jalpha18 Purified**
Cat. Number/Size **11-927-C100 0.1 mg**
For Research Use Only.
Not for use in diagnostic or therapeutic procedures.

Antigen	TCR Valpha24-Jalpha18
Clone	6B11
Format	Purified
Reactivity	Human, Non-human primates
Application	ICC, FC (QC tested), IHC(F)
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody 6B11 recognizes an extracellular epitope within V alpha 24 J alpha 18 TCR, which is expressed by invariant NK T cells (iNKT).

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.