

Technical Data Sheet

Product

**Anti-HRP Purified**

Cat. Number/Size

**11-262-C025**

**0.025 mg**

**11-262-C100**

**0.1 mg**

**For Research Use Only.**

**Not for use in diagnostic or therapeutic procedures.**

Antigen

HRP

Clone

HP-03

Format

Purified

Reactivity

Horseradish, Tagged fusion proteins in all species

Application

WB, ICC, ELISA

Application details

Immunocytochemistry: The antibody HP-03 was successfully used for staining of formaldehyde-fixed, Triton-permeabilized cells transfected with HRP gene.

Isotype

Mouse IgG1

Specificity

The antibody HP-03 reacts with horseradish peroxidase. It is suitable for PAP detection system.

Immunogen

Horseradish peroxidase (Grade I)

Concentration

1 mg/ml

Preparation

Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).

Formulation

Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage and handling

Store at 2-8°C. Do not freeze.

Images and References

[www.exbio.cz](http://www.exbio.cz)

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.