

Technical Data Sheet

Product	Anti-NTAL Purified	
Cat. Number/Size	11-374-C025	0.025 mg
	11-374-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	NTAL
Clone	NAP-07
Format	Purified
Reactivity	Human, Mouse
Application	FC, IP, WB, ICC
Isotype	Mouse IgG1
Specificity	The antibody NAP-07 reacts with an intracellular epitope of Non-T cell activation linker (NTAL), also known as LAB (linker of activated B cells), a 25 - 30 kDa transmembrane adaptor protein associating with lipid raft components in the plasma membrane of B lymphocytes, NK cells and myeloid cells.
Other names	LAB, LAT2
Immunogen	Recombinant cytoplasmic domain (aa 91-243) of human NTAL.
Entrez Gene ID	56743
Gene name	Lat2
NCBI Full Gene Name	linker for activation of T cells family member 2
UniProt ID	Q9JHL0

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

Images and References	www.exbio.cz
-----------------------	--

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.