



	Technical Data Sheet	
Product	Anti-Hu CD155 Biotin	
Cat. Number/Size	1B-851-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD155	
Clone	SKII.4	
Format	Biotin	
Reactivity	Human	
Application	WB, FC (QC tested), IP	
Application details	Flow cytometry: Recommended dilution: 0.5-4 µg/ml.	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody SKII.4 recognizes an extracellular epitope of CD155, a transmembrane glycoprotein, known as poliovirus and cytomegalovirus receptor.	
Other names	PVR, HVED, NECL5, PVS, TAGE4, CMVR	
Immunogen	human neuroblastoma SK-L cells	
Entrez Gene ID	5817	
Gene name	PVR	
NCBI Full Gene Name	PVR cell adhesion molecule	
UniProt ID	P15151	
Concentration	1 mg/ml	
Preparation	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.