

Technical Data Sheet

Product

Anti-CD49f Purified

Cat. Number/Size

11-913-C100

0.1 mg

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	CD49f
Clone	GoH3
Format	Purified
Reactivity	Non-human primates, Horse, Pig, Rabbit, Cow, Dog, Cat, Human, Sheep, Mouse
Application	FC (QC tested), ICC, IP, IHC-F
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Rat IgG2a
Specificity	The rat monoclonal antibody GoH3 recognizes an extracellular epitope of CD49f, an alpha6 integrin subunit.
Other names	ITGA6, VLA6alpha, platelet gplc, integrin subunit alpha 6
Workshop	HCDM workshop listed
Immunogen	mouse mammary tumor cells
Entrez Gene ID	3655
Gene name	ITGA6
NCBI Full Gene Name	integrin subunit alpha 6
UniProt ID	P23229

Concentration	1 mg/ml
Preparation	Purified by protein-G affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.