



	Technical Data Sheet	
Product	rApi g 1 (Apium graveolens 1.0101)	
Cat. Number/Size	10-P149-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Format	Purified Azide Free	
Application	FC (verified), ELISA	
Description	Recombinant protein rApi g 1 is expressed in Escherichia coli. DNA sequence encoding 167 AAs was fused with His-tag at the C-terminus. A calculated molecular mass of recombinant protein is 17,9 kDa.	
Contents	Recombinant protein rApi g 1 - 1 vial (0.1 mg).	
UniProt ID	P49372	
Concentration	1 mg/ml	
Preparation	Purified by sequential steps of ion exchange and affinity chromatography.	
Formulation	HEPES buffered saline (HBS), pH 7.4	
Storage and handling	Store at -20°C to -80°C.	
Images and References	www.exbio.cz	

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic use. In vivo diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s. © 2019 EXBIO Praha, a.s. All rights reserved.