

Technical Data Sheet

Product

Anti-HLA-DR/DP Purified

Cat. Number/Size

11-240-C025

0.025 mg

11-240-C100

0.1 mg

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	HLA-DR/DP
Clone	HL-40
Format	Purified
Reactivity	Human, Pig
Application	FC (QC tested), IP, IHC(P)
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml.
Isotype	Mouse IgG2a
Specificity	The antibody HL-40 reacts with an extracellular common epitope on beta-chain of human HLA-DR and HLA-DP. DR and DP are the isotypes of human MHC Class II molecules expressed on antigen-presenting cells (APC; dendritic cells, B lymphocytes, monocytes, macrophages).
Immunogen	Raji Burkitt's lymphoma cell line

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.