



	Technical Data Sheet	
Product	Anti-Hu CD49a Purified	
Cat. Number/Size	11-873-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD49a	
Clone	TS2/7	
Format	Purified	
Reactivity	Human, Non-human primates	
Application	FC (QC tested), IP, IHC-F	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody TS2/7 recognizes an extracellular epitope of CD49a, the alpha 1 chain of VLA integrin complex, that serves as an adhesion receptor for collagen IV and laminin 1.	
Other names	alpha 1 integrin, VLA-1 alpha, ITGA1	
Immunogen	human CTL line	
Entrez Gene ID	3672	
Gene name	ITGA1	
NCBI Full Gene Name	integrin subunit alpha 1	
UniProt ID	P56199	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.