



	Technical Data Sheet	
Product	Anti-Hu CD137 Purified	
Cat. Number/Size	11-143-C025	0.025 mg
	11-143-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD137	
Clone	4B4-1	
Format	Purified	
Reactivity	Human, Non-human primates	
Application	FC (QC tested), IP, ELISA	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml	
Isotype	Mouse IgG1 kappa	
Specificity	The mouse monoclonal antibody 4B4-1 recognizes an extracellular conformational epitope on CD137, an approximately 40 kDa type I transmembrane protein of the TNFR family expressed mainly on activated T cells.	
Other names	TNFRSF9, 41BB, 4-1BB, MGC2172	
Workshop	HLDA VI: WS Code C-7	
Immunogen	Recombinant human CD137 ectodomain	
Entrez Gene ID	3604	
Gene name	TNFRSF9	
NCBI Full Gene Name	TNF receptor superfamily member 9	
UniProt ID	Q07011	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.