



	Technical Data Sheet	
Product	Anti-Hu CD106 Purified	
Cat. Number/Size	11-606-C025	0.025 mg
	11-606-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD106	
Clone	STA	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), IP, IHC-F, ELISA	
Application details	Flow cytometry: Recommended dilution: 4-6 µg/ml; positive control: TNF-alpha activated HUVEC cells. Immunohistochemistry (frozen sections): Acetone fixation. ELISA: Capture mAb for soluble CD106.	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody STA recognizes an extracellular epitope of CD106 antigen (VCAM-1), a 100-110 kDa type I membrane protein of the immunoglobulin superfamily, a crucial mediator of leukocyte adhesion, and a costimulation molecule.	
Other names	VCAM-1, INCAM-100	
Workshop	HLDA V: WS Code A013	
Immunogen	Human DS6 T cell line	
Entrez Gene ID	7412	
Gene name	VCAM1	
NCBI Full Gene Name	vascular cell adhesion molecule 1	
UniProt ID	P19320	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.