

Technical Data Sheet

Product **Anti-Hu TROP2 Purified**
Cat. Number/Size **11-898-C100** **0.1 mg**
For Research Use Only.
Not for use in diagnostic or therapeutic procedures.

Antigen	TROP2
Clone	TrMab-6
Format	Purified
Reactivity	Pig, Human
Negative species	Mouse
Application	WB, FC (QC tested), IHC, ICC
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml. Western blotting: Recommended dilution: 1 µg/ml. Immunocytochemistry: Recommended dilution: 10 µg/ml; the cells can be fixed with 4% PFA and permeabilized with 0.1% Triton X-100 before antibody staining, when needed. Immunohistochemistry (paraffin sections): Recommended dilution: 10 µg/ml. Heat-mediated antigen retrieval in citrate buffer, pH = 6.
Isotype	Mouse IgG2b kappa
Specificity	The mouse monoclonal antibody TrMab-6 recognizes an extracellular epitope of TROP2, a type I transmembrane glycoprotein. This antibody is useful for detection of TROP2 in breast cancer.
Other names	TACSTD2, EGP1, GP50, M1S1, GA733-1
Immunogen	TROP2-transfected CHO cells
Entrez Gene ID	4070
Gene name	TACSTD2
NCBI Full Gene Name	tumor associated calcium signal transducer 2
UniProt ID	P09758

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.