

Technical Data Sheet

Product	<b>Anti-TRIM Purified</b>	
Cat. Number/Size	<b>11-344-C025</b>	<b>0.025 mg</b>
	<b>11-344-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	TRIM
Clone	TRIM-04
Format	Purified
Reactivity	Human
Negative species	Mouse
Application	FC, IP, WB
Isotype	Mouse IgG2a
Specificity	The antibody TRIM-04 recognizes an intracellular epitope of T cell receptor-interacting molecule (TRIM), a 30 kDa adaptor protein expressed by T cells.
Other names	TRAT, TCRIM
Immunogen	Recombinant intracellular domain (aa 29-186) of human TRIM.
Entrez Gene ID	50852
Gene name	TRAT1
NCBI Full Gene Name	T cell receptor associated transmembrane adaptor 1
UniProt ID	Q6PIZ9

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.