

Technical Data Sheet

Product	<b>Anti-Hu CD172g Purified</b>	
Cat. Number/Size	<b>11-881-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	CD172g
Clone	OX-119
Format	Purified
Reactivity	Human
Application	FC (QC tested), IP
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody OX-119 recognizes an extracellular epitope on CD172g, an approximately 55 kDa transmembrane glycoprotein expressed on most T cells, as well as on NK cells and some B cell populations.
Other names	SIRP gamma, SIRP-B2
Immunogen	recombinant human CD172g
Entrez Gene ID	55423
Gene name	SIRPG
NCBI Full Gene Name	signal regulatory protein gamma
UniProt ID	Q9P1W8

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.