



	Technical Data Sheet	
Product	Anti-Hu CD3 Purified	
Cat. Number/Size	11-381-C025	0.025 mg
	11-381-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD3	
Clone	MEM-92	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), IP	
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml	
Isotype	Mouse IgM	
Specificity	The antibody MEM-92 reacts with an extracellular epitope on epsilon chain of human CD3 complex, a part of a bigger multisubunit complex of the T cell receptor (CD3/TCR) expressed on peripheral blood T lymphocytes and mature thymocytes.	
Other names	CD3E, T3E, TCRE	
Workshop	HLDA IV: WS Code T 97	
Immunogen	Human peripheral blood lymphocytes.	
Entrez Gene ID	916	
Gene name	CD3E	
NCBI Full Gene Name	CD3 antigen, epsilon polypeptide	
UniProt ID	P07766	
Concentration	1 mg/ml	
Preparation	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).	
Formulation	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.