

Technical Data Sheet

Product	<b>Anti-Ag85b (M. tuberculosis) Purified</b>	
Cat. Number/Size	<b>11-189-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	Ag85b (M. tuberculosis)
Clone	PAb (189)
Format	Purified
Reactivity	Mycobacterium tuberculosis
Application	WB, ELISA
Isotype	Rabbit polyclonal
Specificity	The rabbit polyclonal antibody to Ag85b reacts with the Mycobacterium tuberculosis protein Ag85b (FbpB), which is a diacylglycerol acyltransferase/mycolyltransferase.
Other names	FbpB, mycolyl transferase 85B
Immunogen	Recombinant Mycobacterium tuberculosis Ag85b produced in E. coli
Entrez Gene ID	885785
Gene name	fbpB
NCBI Full Gene Name	diacylglycerol acyltransferase/mycolyltransferase
UniProt ID	I6YBV4

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.