

Technical Data Sheet

Product **Anti-Hu IL-17A Purified**  
Cat. Number/Size **11-937-C100** **0.1 mg**  
**For Research Use Only.**  
**Not for use in diagnostic or therapeutic procedures.**

---

Antigen	IL-17A
Clone	9F9
Format	Purified
Reactivity	Human
Application	IP, ELISA, FC (QC tested), ICC
Application details	Flow cytometry: Recommended dilution: 0.5-4 µg/ml. Intracellular staining.
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody 9F9 recognizes human interleukin 17A (IL-17A; secreted or intracellular).
Other names	Interleukin 17
Immunogen	mammalian-derived human IL-17-IgG fusion protein, boost with recombinant human IL-17A
Entrez Gene ID	3605
Gene name	IL17A
NCBI Full Gene Name	interleukin 17A
UniProt ID	Q16552

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.