



Technical Data Sheet

Product	Anti-CPNE7 Purified	
Cat. Number/Size	11-542-C025	0.025 mg
	11-542-C100	0.1 mg
For Research Use Only.		
Not for use in diagnostic or therapeutic procedures.		
Antigen	CPNE7	
Clone	CPNE7-01	
Format	Purified	
Reactivity	Human	
Application	WB	
Isotype	Mouse IgG	
Specificity	The mouse monoclonal antibody CPNE7-01 recognizes C terminus of human CPNE7, a calcium-binding cytoplasmic protein expressed mainly in brain.	
Other names	copine 7	
Immunogen	Bacterially produced GST-fused C terminus of human CPNE7.	
Entrez Gene ID	27132	
Gene name	CPNE7	
NCBI Full Gene Name	copine 7	
UniProt ID	Q9UBL6	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.