



	Technical Data Sheet	
Product	<b>Anti-CTBP1 Purified</b>	
Cat. Number/Size	<b>11-511-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	
Antigen	CTBP1	
Clone	PAb (511)	
Format	Purified	
Reactivity	Mouse, Human	
Application	IP, WB, ICC	
Application details	Western blotting: Reducing conditions. Immunocytochemistry: Recommended dilution 1 µg/ml, acetone / methanol fixation.	
Isotype	Rabbit polyclonal	
Specificity	The polyclonal antibody reacts with CtBP1, a 49 kDa transcriptional corepressor (intracellular antigen), that plays a critical role in regulating gene expression.	
Other names	BARS, HADDTS	
Immunogen	Bacterially expressed full-length human CtBP1	
Entrez Gene ID	1487	
Gene name	CTBP1	
NCBI Full Gene Name	C-terminal binding protein 1	
UniProt ID	Q13363	
Concentration	1 mg/ml	
Preparation	Purified by ligand affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.