



	Technical Data Sheet	
Product	Anti-PRKAR2A Purified	
Cat. Number/Size	11-428-C025	0.025 mg
	11-428-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	PRKAR2A	
Clone	Hs-36	
Format	Purified	
Reactivity	Human	
Application	ICC, WB	
Application details	Immunocytochemistry: Recommended dilution: 10 µg/ml; membrane permeabilization (acetone) is essential.	
Isotype	Mouse IgM	
Specificity	The antibody Hs-36 reacts with PRKAR2A (protein kinase A regulatory type II alpha subunit), an intra-acrosomal protein.	
Other names	PKAR2A, PKR2, PKAR2, proteinkinase A regulatory 2	
Immunogen	Freshly ejaculated human sperms were washed in PBS and extracted in 3% acetic acid, 10% glycerol, 30 mM benzaminidine. The acid extract was dialyzed against 0.2% acetic acid and subsequently used for immunization.	
Entrez Gene ID	5576	
Gene name	PRKAR2A	
NCBI Full Gene Name	protein kinase cAMP-dependent type II regulatory s	
UniProt ID	P13861	
Concentration	1 mg/ml	
Preparation	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).	
Formulation	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.