



	Technical Data Sheet	
Product	Anti-Hu CD253 Purified	
Cat. Number/Size	11-316-C025	0.025 mg
	11-316-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD253	
Clone	2E5	
Format	Purified	
Reactivity	Human	
Negative species	Mouse	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml.	
Isotype	Mouse IgG1	
Specificity	The antibody 2E5 reacts with an extracellular epitope within C-terminal half of TRAIL (APO-2L), a 21 kDa cytotoxic protein, activator of rapid apoptosis in tumor cells. TRAIL is mainly expressed in spleen, lung, prostate and also in many other tissues.	
Other names	TRAIL, Apo-2 ligand, Apo-2L, TNFSF10, APO2L, TNLG6A	
Workshop	HLDA IX	
Immunogen	Recombinant soluble fragment (aa 95-281) of human TRAIL.	
Entrez Gene ID	8743	
Gene name	TNFSF10	
NCBI Full Gene Name	TNF superfamily member 10	
UniProt ID	P50591	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.